

# Microscope Body Specifications

Available Zoom Bodies Range:	0.8X - 5X, 0.8X - 6.4X, 0.8X - 8X
Zoom Ratio:	1:6, 1:8, 1:10 (respectively)
Magnification w/ 10X eyepieces	8X - 50X, 8X - 64X, 8X - 80X (respectively)
Field of View	28.8mm-4.6mm, 28.8mm-3.6mm, 28.8mm-2.9mm (respectively)
Working Distance	78mm
Eyepieces Inclination:	20° and 0° - 45°
Dioptric Adjustment:	±5 Dual Independent
Interpupillary Adjustment:	51mm - 75mm
Microscope Ring Mount Diameter:	76mm
Available Eyepieces:	10X, 20X and 30X
Available Objective Lens:	0.5X Achromatic, 0.5X Plan Apo, 1X Plan, and 2X Plan
Available Video Coupler:	0.35X and 0.5X
Available Video Slitter:	Single bean port
Optional Iris:	Yes